Instructions for use

General

Henry Schein® Eurosil® Condensation Cured Silicone Correction Material (Henry Schein® Eurosil® Correction Material) is a low viscosity impression material based on condensation curing silicones. Henry Schein® Eurosil® Correction Material shows hydrophilic properties, thus improving wettability and leading to better reproduction of subgingival areas.

Classification

Henry Schein® Eurosil® Correction Material meets the requirements of: DIN 13913 2-c  
ISO 4823 Cat. B, Type 3 (low viscosity)

Indication

Correction impression material for the double-impression technique in conjunction with Henry Schein® Eurosil® Putty.

Dispensing and mixing

Place equal lengths (10:1 by weight or 11.7:1 by volume) of Henry Schein® Eurosil® Correction Material and Henry Schein® Eurosil® Catalyst on the mixing pad and blend with a spatula. Mix at least 30 seconds. Color must be homogeneous.

Impression taking

The ready mixed material is loaded into the tray. A part of the mixture might be additionally syringed directly onto the teeth. Insert the tray not later than 2 minutes after start of mix into the mouth. Maintain position for at least 3½ minutes.

Note: Times mentioned are valid for 23°C (73°F). Higher temperature shortens, lower temperature prolongs these times.

Casting

The impression should be stored under normal room conditions (max. 25°C/77°F) and poured within 2-4 hours.

Electroplating

All impressions can be copper- or silver-plated.

Disinfection

Impressions can be disinfected with standard disinfection solution.

Storage

The shelf life is at least 3 years, when stored at 23°C (73°F). Do not store above 25°C (77°F).

Additional notes

Avoid contact of Henry Schein® Eurosil® Correction Material with eyes. If contact occurs, wash the eyes with copious amounts of water and seek ophthalmological care.

Ready cured Henry Schein® Eurosil® Correction Material is resistant against solvents. Avoid spots on laboratory coats and other clothing.

Technical data

Mixing time .....	30 sec.
Total working time (incl. mixing time).....	2 min.
Setting time (incl. mixing time).....	4 min.
Recommended time in mouth.....	3½ min.
Dimensional change after 24 hours.....	less than 0.8%
Compression set .....	less than 0.3%
Strain in compression .....	10.0%

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Hergestellt in den USA/Gefabriceerd in de U.S.A.